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For dental use only

RECIPROC® instruments for root canal preparation

1. Composition

The working part of the RECIPROC® instruments is made of M-Wire® nickel-titanium. The M-Wire® alloy undergoes a thermal-treatment process for more flexibility and a greater cyclic fatigue resistance, compared to traditional nickel-titanium.

2. Contents

RECIPROC® blister package, 4 or 6 sterile instruments each

3. Indications for Use

- Preparation of the root canal system.

RECIPROC® is used in reciprocation with a corresponding motor with specific RECIPROC® motor settings, e.g. VDW.SILVER® RECIPROC®. Due to its specific design, the instrument cuts counterclockwise in reciprocation.

4. Contraindications

In some canals, the hand file used for determining the working length (after the RECIPROC® instrument has reached 2/3 of the estimated working length) may need to be pre-curved in order to reach the working length. This indicates the presence of an abrupt apical curvature. The use of RECIPROC® instruments is contra-indicated in this instance. In these cases the canal preparation has to be finished with hand files. This limitation also applies to continuous rotary instruments.

5. Warnings

None

6. Precautions

- The instruments may only be used by qualified persons in dentists' offices or clinics.
- A RECIPROC® instrument is designed for single use in maximum one molar and cannot be reused. The instrument cannot be sterilised.
- Examine the instrument during treatment for signs of wear such as untwisting, bending or damaged cutting edges. In these cases, instruments should be discarded.

- Clean the debris from the flutes, every three pecks, in an interim stand. Irrigate the canal according to the appropriate irrigation protocol.
- Only light pressure should be applied on RECIPROC® instruments. If resistance is encountered or if the instrument does not advance, clean the instrument, irrigate the canal according to the appropriate irrigation protocol and check the canal for possible blockage with a C-PILOT® file ISO size 10. Continue with the treatment.
- Use a rubber dam.

7. Adverse Reactions

None

8. Step-by-Step Instructions - Preparation with RECIPROC®

- **Create a straight line access to the root canal entrance**
- **Select the correct RECIPROC® instrument:**
In most cases, the R25 will be suitable in size for the root canal treatment. Consult the pre-operative radiograph to see if the canal is likely to be considered as narrow, medium or wide.

If the canal is partially or completely invisible on the radiograph:
The canal is considered **narrow = use an R25**.

If the canal is completely visible on the radiograph:

1. Take an ISO size 30 hand instrument; insert it passively into the canal. If it reaches working length, the canal is considered **large = use the R50**.
2. If an ISO size 30 hand instrument does **not** go passively to working length, try an ISO size 20 hand instrument. If this goes passively to length, the canal is considered **medium = use an R40**.
3. If a size 20 hand instrument does **not** go passively to length, the canal is considered **narrow = use an R25**.

Passively means that the instrument goes directly to working length with a gentle watch winding movement (small right left rotations) but without filing action.

- **Preparation with the R25 – Getting started:**

Prior to preparation with an R25, estimate the working length with the help of a pre-operative radiograph. Set the stopper at 2/3 of that length.

- **Preparation with the R40 and R50 – Getting started:**

Before selecting the appropriate RECIPROC® instrument, check the working length by using a C-PILOT® or K-file. Set the stopper at that length.

- **Further steps:**

1. Place irrigant in the access cavity of the root canal.
2. Introduce the RECIPROC® instrument into the canal. Press the motor foot pedal when the instrument is at the root canal orifice.
3. Move the instrument in a slow in-and-out pecking motion. The amplitude of the in-and-out movements should not exceed 3 mm. Only very light pressure should be applied. The instrument will advance easily in the canal.
One in-and-out movement = **one peck**.
If necessary, the instrument can be used in a brushing motion to enlarge the root canal entrance.
4. **Clean the instrument in the interim stand after 3 pecks.**
5. **Irrigate the canal.**
6. **Make sure the canal is free with an ISO size 10 C-PILOT® file.**
7. Repeat step 3-6 until approx. 2/3 of the working length has been reached (indicated by the silicone stopper).

8. **When using an R25:**

Electronic length determination of the working length.

When using a R40 or R50:

After preparation of 2/3 of the working length, the working length should be rechecked with an apex locator.

9. Afterwards continue with step 3-6 until full working length has been reached.

Further languages are available upon request.

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VDW GmbH