HELIX TEST PROCESS CHALLENGE DEVICE

INSTRUCTIONS

1. Open the Teflon ® tube
2. Fold the strip at the incision. The ink indicator must be turned inwards.
3. Insert the strip thus folded into the cannula. The folded side should be toward the inside of the chamber. The test results might not be distorted by the incorrect positioning of the strip.
4. Close the Teflon ® tube.
5. Make a sterilization cycle. The parameters necessary for sterilization to obtain an ink color change indicator are:
134 ° C saturated steam for 3.5 minutes
121 ° C saturated steam for 15 minutes
6. When the test open the capsule of the device and pull the indicator.
7. Where all the air inside the device has been effectively removed, the saturated steam has penetrated inside the tube until it reaches the chamber and the temperature of 134 ° C was maintained for at least 3.5 minutes (in Alternative: 121 ° C for 15 minutes) the ink indicator will be tacked from yellow to black, indicating that all parameters of sterilization have been met.

The HELIX TEST SYSTEM complies with standards: EN 867-5 / ISO 11140-1 Class 6

HELIX TEST PROCESS CHALLENGE DEVICE

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